
# Product Innovations 2023

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All rights reserved.

Reproduction, also by extract and reproducing photos are only permitted with written authorization of EDENTA AG, Switzerland / Edenta Etablissement, Liechtenstein. All offers, orders and deliveries are subject to Edenta's „General sales and Delivery Terms“. We reserve the right to modify our range of products and their design as well as to deviate from the illustrations and data shown.

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<b>Zahnarzt-Praxis / Dental Practice / Cabinet Dentaire</b>	<b>Page</b>
Bone Cutter C162AAC.316.016 - C162AAC.317.016	4
Bone Cutter C162SAC.316.016	5
Bone Cutter C254AAC.316.012	6
Surgery Diamond FGXXL - G254T.317.014	7
Surgery Diamond FGXXL - D254T.317.016	8
Osteotom FGL - GSD1T.315.016 - GSD1T.315.018 - GSD1T.315.021	9
Crown Cutter C4AKXL.315.014 - C4AKXXL.315.014	10
Crown Cutter C34LT.315.012 - C34LT.315.016	11
KT-Diamond Turbo FG - Crown Cutter - KT856.314.016 - KT881.314.016	12
Adhesive Remover T379X.204.018 - T379X.204.023	13
Diamonds WhiteTIGER with gold plated Shank (FG)	14 - 15
Diamond Strips KFO	16
StarTec Lamellae Polisher RA for polishing Lithium Disilicate (LS2)	17
 <b>Dental-Labor / Dental Laboratory / Laboratoire Dentaire</b>	
CeraGloss Lamellae Polisher HP for polishing All-Ceramics	18
StarTec Lamellae Polisher HP for polishing Lithium Disilicate (LS2)	19
StarGloss Lamellae Polisher HP for polishing Zirconia	20
CeraTec Wheel for preparing High-Performance Zirconia - 953.130HP	21
Tungsten Carbide Cutter toothing VZ26 - 7226.060HP	22
Tungsten Carbide Cutter toothing VZ27 - 7727.029HP	23
Cutters for CAD/CAM systems Amann Girrbach	24
Cutters for CAD/CAM systems Roland	25 - 26
Cutters for CAD/CAM systems VHF	27

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## Bone cutter C162AAC

Conical bone cutter used for preparing bone and tooth structure. The Lindemann blades, newly optimised with ZrN (zircon-nitride) hard coating, guarantee quiet instrument running at the highest cutting capacity and perfect cutting performance for enhanced cutting sensation and maximum cutting control. The hard ZrN (zircon-nitride) coating ensures the cutter meets the high demands placed on a high-quality surface during surgical bone preparation.

**Application:** Separating teeth and tooth roots, exposing and removing impacted wisdom teeth, removing of root remains.

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
## Bone cutter C162SAC

Bone cutter for gentle trimming of bone and tooth structure. The special S-cut for optimum cutting sharpness and perfect blade action for an even better cutting sensation and maximum cutting control.

The hard ZrN (zircon-nitride) coating ensures the cutter meets the high demands placed on a high-quality surface during surgical bone preparation.

Indications: Exposing and removing retained teeth, removing of root remains, preparation of bone structure and bone lids, apicectomies, cutting teeth and tooth roots.

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## Bone cutter C254AAC

Bone cutter used for minimally invasive preparation of bone and tooth structure. The special A-cut blades guarantee very quiet instrument running at the highest cutting capacity and maximum conservation of bone structure to be prepared. The long, slimline instrument neck provides excellent visibility during preparation. The hard ZrN (zircon-nitride) coating ensures the cutter meets the high demands placed on a high-quality surface during surgical bone preparation.

**Application:** Exposing and removing impacted wisdom teeth, apicectomies, preparation of bone structure and bone lids, separating teeth and tooth roots.

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
## Surgical diamond G254T

Surgical diamond instrument for removing and smoothing bone surfaces as well as for separating teeth and effectively cutting tooth substance in dental osteotomy. The multilayer diamond coating ensures gentle treatment of the bone structure and therefore also less heat development. In addition, the high-performance TiN coating ensures optimum removal of bone and tooth substance as well as a long service life of the diamond instruments.

**Grit size:** blue / standard = 126 µm – ISO 524. Ensure sufficient cooling during the work process.

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
## Surgical diamond D254T

Surgical diamond instrument for atraumatically separating teeth and effectively cutting tooth substance in root tip resection. The multilayer diamond coating ensures gentle treatment of the tooth substance and therefore also less heat development. In addition, the high-performance TiN coating ensures optimum removal of tooth substance as well as a long service life of the diamond instruments.

Not recommended for chipping bone.

**Grit size:** green / coarse = 151 µm – ISO 534. Ensure sufficient cooling during the work process.

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## Osteotome GSD1T

For removing and smoothing bone as well as cutting tooth roots. The special shape together with the multilayer diamond coating ensure a gentle fit to the bone structure and thus less heat development. In addition, the high-performance TiN coating ensures optimum removal of bone and tooth substance as well as a long service life.

**Grit size:** green/coarse = 151 µm – ISO 534. Ensure sufficient cooling during the work process.

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## Crown cutter C4AK

For cutting low-fusing veneering ceramics and metal frameworks in a single working step, easily and quickly. The novel MultiCut toothing enables fast chip removal and ensures good instrument stability with a high breaking strength. The new toothing system ensures smooth and low-vibration operation. The crown cutters are to be placed at an angle of 45°.

Caution: the crown cutters must not be used for separating all-ceramic restorations made of extremely hard ceramics, e.g. zirconium oxide. Use the K-diamonds for separating all-ceramic materials, as they have been developed specifically for this application. Speed ideally 160,000 rpm, contact pressure: 0.5 N

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

## Crown cutter C34LT

New generation of crown cutters with TiN high-performance coating for optimal chip removal and smooth, low-vibration operation. The toothing with crosscut and overhead cutting edge ensures maximum cutting performance as well as a long service life, even when cutting non-precious metal crowns. The crown cutter C34LT.315.012 is to be placed at an angle of 45°. The crown cutter C34LT.315.016 can be used both at a 45° angle as well as flat.

Notes on use: 1. cut open the ceramic veneer with a diamond instrument (K-diamond) 2. cut the metal framework with the crown cutter – optimal speed 160,000 rpm, contact pressure 0.5 N

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## Crown cutter KT856 / KT 881

For cutting all-ceramic restorations such as zirconium oxide and other all-ceramic materials. The exceptional stability of the special diamond grit provides for a long service life of the KT diamond instruments. The multilayer electroplated bond prevents the diamond grit from loosening prematurely. This results in an optimal removal rate with low vibrations when cutting all-ceramic materials.

The crown cutter KT856 is to be placed at an angle of 45°. To remove an all-ceramic restoration, the restoration must be slit along its axial wall towards the centre of the occlusal surface or the incisal edge. The restoration is then bent open with a suitable instrument (lever), causing it to fracture. The KT881 is suitable for grinding off fragment residues. Speed ideally 160,000 rpm with the red contra-angle, contact pressure: < 2 N

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## Adhesive remover T379X

Tungsten carbide instruments for the removal of adhesive residues after removal of orthodontic brackets. Used at low speed, the instruments are highly effective for the coarse removal of remaining adhesive and only create a slight roughness depth. The special cutting geometry enables targeted and low-vibration removal of soft adhesive residues. The egg shape is particularly suitable for palatal removal of the adhesive when applying the lingual technique.

Advantages of the adhesive remover: a special TiN coating prevents smearing of the instruments' cutting edges. This allows for smooth and controlled removal of the adhesive residues, with only minimal heat development. Speed ideally 30,000 rpm.

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

## WhiteTiger diamonds

New with gold-plated shaft. This makes it very easy to distinguish the high-quality WhiteTiger diamond instruments from conventional instruments. A special bond securely binds the multilayer diamond particles. The additional, wear-resistant LLD coating ensures an improved service life as well as increased abrasion performance compared to conventional diamond instruments. WhiteTiger diamonds thus enable gentle treatment. The multilayer diamond structure ensures fast and effective material abrasion. The LLD coating reduces clogging of the chip spaces and thus lowers heat generation. Specifically for crown abutment and cavity preparation. The increased abrasion rate and better cooling compared to conventional instruments deliver an efficient and at the same time gentle, optimal preparation result.

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














## Application:

### Crown abutment and cavity preparation.

The multilayer diamond coating together with the LLD coating ensure increased material removal with good chip flow for cooler grinding. Compared to conventional diamond instruments, this provides higher tooth substance removal to deliver an efficient and gentle preparation result.

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## Orthodontic diamond finishing strips

For interproximal enamel reduction (IPR), also referred to as “stripping“ or “recontouring“. IPR is a technique in which part of the enamel layer is removed to create sufficient space for orthodontic alignment.

Diamond-coated on both sides, the diamond finishing strips enable simple and precise reduction of the interproximal tooth substance. The flexible strips follow the natural contours of the teeth to provide optimal surface design.

The diamond finishing strips are available in 4 different thicknesses: 0,1 / 0,2 / 0,3 und 0,4 mm. Three different grit sizes are available, according to the colour coding; yellow for superfine grit, red for fine grit and blue for medium grit. The stainless steel strips can be autoclaved and reused.

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## StarTec RA Porcelain Polishers

2-Step Polishing system with diamond particles, used for producing an intraoral high-lustre polish on zirconium dioxide, lithium disilicate (LS2), zircon-reinforced lithium silicate (ZLS) as well as hybrid-ceramic. The highly flexible polishing lamellae adapt to any surface structure thus creating perfect polishing results, even on areas difficult to polish such as occlusal surfaces. Recommended speed 6.000 rpm. Always use coolant, 50L/min, when polishing. External cooling should be used when polishing with lamellae polishers, which should be directed between the polisher and ceramic filling.

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

## CeraGloss HP polisher for ceramics

Two-step lamella polishing system with diamond grit, for high-lustre polishing of all-ceramic/ and ceramic materials. The individual lamellae support each other during the polishing process while retaining their flexibility. Therefore, the polishing lamellae adapt to every surface structure and create perfect polishing results. For high-lustre polishing of occlusal, labial, incisal, buccal and interproximal areas.

Polishing step 1 (blue) for surface smoothing and preparation for final polishing.  
 Polishing step 2 (yellow) for high-lustre polishing without additional glaze firing  
 Recommended speed 5.000 rpm.

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## StarTec HP polisher for ceramics

Two-step lamella polishing system with diamond grit, for high-lustre polishing of zirconium dioxide, lithium disilicate, zirconium reinforced lithium silicate (ZLS), hybrid ceramics, as well as all common veneering ceramics. The individual lamellae support each other during the polishing process while retaining their flexibility. Therefore, the polishing lamellae adapt to every surface structure and create perfect polishing results. For high-lustre polishing of occlusal, labial, incisal, buccal and interproximal areas.

Polishing step 1 (purple) for surface smoothing and preparation for final polishing.  
 Polishing step 2 (yellow) for high-lustre polishing without additional glaze firing  
 Recommended speed 5.000 rpm.

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## StarGloss HP polisher for ceramics

Two-step lamella polishing system with diamond grit, for high-lustre polishing of ceramic and all-ceramic materials such as zirconium oxide ZrO<sub>2</sub>, aluminium oxide. The individual lamellae support each other during the polishing process while retaining their flexibility. Therefore, the polishing lamellae adapt to every surface structure and create perfect polishing results. For high-lustre polishing of occlusal, labial, incisal, buccal and interproximal areas.

Polishing step 1 (pink) for surface smoothing and preparation for final polishing.  
 Polishing step 2 (grey) for high-lustre polishing without additional glaze firing  
 Recommended speed 5,000 rpm.

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### CeraTec ceramic grinder

Abrasive trimmer with ceramic bond and specially matched diamond grit for processing zirconium materials without heating the all-ceramic material structure. For the efficient and material-conserving processing of high-performance zirconium, lithium disilicate, zirconium reinforced lithium silicate (ZLS), hybrid ceramics, resin nano ceramics (RNC) as well as all common veneering ceramics. CeraTec wheel (953.130HP), ideal for smoothing attachment points / retaining pin attachments and larger surfaces. The small size of the wheel also has a positive effect on the field of vision. Optimal speed 5,000 rpm.

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## Economical machining of titanium and titanium alloys with coarse titanium toothing VZ26

Tungsten carbide milling cutter with VZ26 toothing for machining titanium and titanium alloys. The high toughness as well as the low modulus of elasticity lead to increased heat generation and increased wear of the machining tools when cutting titanium materials. The special coarse titanium toothing VZ26 with its low number of cutting edges results in high removal rates as well as a prolonged service life of the milling cutter. The deep crosscut ensures easy penetration into the tough titanium materials. The large chip spaces prevent clogging of the cutting edges. For use in the crown-bridge technique.

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## Economic preparation of non-precious metal alloys with fine titanium blades VZ27

Generally, more material has to be removed during preparation of non-precious metals in the dental laboratory because of technical reasons during casting. The preparation instruments are therefore subjected to a large amount of stress and have only a short service life. New nonprecious metal alloys are slightly softer, but are still difficult to machine. Cutters are consequently required that are as wear resistant as possible and at the same time also ensure high material removal rate. The newly developed fine titanium blade VZ27 has a very high cutting efficiency ensuring a high material removal rate while simultaneously smoothing the surface of the material, which then allows it to be easily polished.

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Certain products and names are brand-, patent-, and copyright-protected company and brand names.