

GLOBAL

SURGICAL™ CORPORATION



*Shown with 1080 HD Video,
Beamsplitter, Carr Binocular
Extender, and Floorstand
Monitor Mount


A-Series™ Microscope Accessories

Global Microscopes Accessories



Multi-Focal Lens

Provides extended focal range of 150mm (6"). User can focus anywhere from 200-350mm (8-14"), all without adjusting the microscope head or the patient. Minimizes learning curve and increases productivity.



Carr Binocular Extender

A 45° wedge that adds additional height and 2.5" (63.5mm) of reach from the microscope head to the user. Helps operator maintain an upright, ergonomic position – alleviating fatigue due to neck and back strain. Great for taller and more robust dentists.



Counterbalance Coupler Arm

Enables accessories to be added while ensuring their weight is counter-balanced. This maintains easy, balanced movement of the microscope head when used with heavier cameras, observation systems, and more.



McDonald Arm 17" Support Extension (43.2cm)

Provides additional reach and maneuverability for larger operatories. Can be used to mount scope above cabinets and for other complex installations.



50/50 and Virtual Beamsplitters

Uses prisms to reflect light from image to accessories such as video or digital cameras. Accessory can be added to one or both sides of the microscope using Single-Port or Dual-Port options.



1080 HD Video Camera*

Provides live feed that can be displayed on a monitor to enhance patient education, assistant's teamwork and to replace an intraoral camera. Can be paired with optional capture device to record video and still images.

*HD camera shown with Beamsplitter and Video Adapter

GLOBAL
SURGICAL™ CORPORATION

www.globalsurgical.com

GSC 17005